



ONEWAY
BIOMED



Two piece implants
DENTAL IMPLANT SYSTEM
GIO, GCO AND GIW

SYMBOLS FOR IMPLANT PROPERTIES AND PROSTHETIC SOLUTIONS

Small Abutment Head



Multi-Unit Abutment



Apical Wide Thread



Large Abutment Head



Micro Thread



For Screw-in Prosthetics



Ball Abutment



Bendable



For Cemented Prosthetics



Localicer® Abutment



ONEWAY
BIOMEED

THE ADVANTAGES

OF THE ENDOSSEOUS DENTAL IMPLANT SYSTEMS **GIO**, **GCO** and **GIW**

GIO, **GCO** and **GIW** implants consist of highly break-proof titanium alloy Ti6Al4V. **GIO** is a successful implant with octagon socket and 8° inner cone for rotationally-secured suprastructures. The surface is roughened in the endosseous region. **GCO** is a precision conical implant with octagon socket and 8° inner cone. Compression screw thread for stable endosseous anchorage. For rotationally-secured suprastructures. **GIW** implants with internal octagon and 8° inner cone are for anti-rotation super-structures.

The prescribed or recommended tightening torques for implants, abutments and screws can be found on our website:

www.implant.com/en/downloads




Secure rotation protection by precision inner octagon


Tightness due to 8° inner cone

Universal application for permanent and removable prosthodontics


Made of highly resistant titanium alloy

Smart instrument tray





Internal octagon



Wide range of sizes:

GIO

Enossal length
Enossal Ø

7 - 17 mm
3.3 - 4.8 mm



GCO

Enossal length
Enossal Ø

7 - 15 mm
3.7 - 4.8 mm



GIW


Enossal length
Enossal Ø

9 - 15 mm
4.8 mm




SURGERY

1.1a. Recommended drill sequences for GIO implants




1.1b. Recommended drill sequences for GCO implants



1.2 Tap for lower jaw (only for use with GIO)

- TAP 1 Size 3.3 (for implants Ø 3.3 mm)
 Size 4.1 (for implants Ø 4.1 mm)
 Size 4.8 (for implants Ø 4.8 mm)



2. Implant packaging GIO/GCO



Original packaging




Open package using the flap. Remove the label and stick it into your patient's record.



The opened package contains the implant in a sterile tube (primary package).

3. Removing the implant from the sterile tube


1. Open the lid. The implant is connected to the lid through a breakable section.
2. Remove the implant without touching the inner walls of the tube.




4. Handling

Assembling the placement aid

Attach the insertion tool to the implant (GIO or GCO) by holding the Lid, to which the implant is secured, with your other hand.



After you have attached the insertion tool, firmly hold the lid in your hand and break the implant off the lid. Insert the implant manually until it fits tightly into the jaw.



5. Insertion

Using the ratchet or contra-angle handpiece: screw the implant clockwise into the bone cavity. The endosseous part of the implant must be **completely** covered by the bone. The polished implant neck is located partially outside of the bone or at bone level. We recommend to screw min. 1 mm of the polished implant neck inside of the bone.




6. Remove insertion tool from implant

Release of the insertion tool or the contra-angle handpiece from the implant:

Using HT 1.25 loosen the screw in the insertion tool while loading the insertion tool IT GCO in insertion direction with the ratchet.

Pull insertion tool off the implant and separate from the contra-angle hand piece.

Attach IT ITV Ratchet adapter with RAT2 ratchet to the ITV insertion tool. Using the HAS flat wrench firmly hold the lower hexagon of the ITV. Remove the ITV from the implant using the ratchet ("Out" position).




6. Result

The implant is fully inserted and ready for healing.




8. Post-operative treatment

Close the implant with a cover screw CST GTI (tool: HT 1.25)



After the healing period: Remove cover screw



9. Impression taking

9.1 Impression taking with perforated, individual impression tray

The long pick up screw must be clearly visible, when the impression try is inserted over the impression post.

The impression material must be removed in the area of the screw access prior to hardening.

9.2 Prior to impression taking

Impression taking with an A silicone® such as Safeprint®.

The use of open or closed impression trys is possible.

9.3 Impression taking





Remove OLT GCO from the implant. OLT GCO remains in the impression.

9.4

View of the impression post OLT GCO in the impression (Pick up method)

9.5


After the impression is taken, the implant is closed with a healing screw (HS) and the impression is sent to the laboratory.



10. Further processing in the laboratory


10.1

Screw IA GCO against the impression post.



10.2

The impression is casted . Then the transfer OLT GCO is unscrewed from the laboratory analogue IA GCO.



10.3

The laboratory analogue is now in the proper position and orientation in the putty.



10.4

Positioning the solid secondary part TLA GCO (straight) or angulated TLA15 GCO, TLA20 GCO or TLA25 GCO, at which the optimal position and adequate angulation must be determined.


Note

The octagon must be completely entered into the analogue. Ensure correct position of the octagon.



10.5

Ensure proper position of the abutment when transferring into the mouth. Tightening torque of the screw during fastening on the implant with: 25 Ncm.



GIO IMPLANTS

Successful implant with octagon socket and 8° inner cone. Connection type: Inner cone + internal octagon / head-Ø: 4.8 mm. Made from titanium alloy ASTM F 136-13/ISO 5832-3. For rotationally-secured supra-structures. Surface is roughened in endosseous region.

		Description	Enossal Ø	Enossal length	REF	Price cat.	
		GIO 3.3 9	3.3 mm	9 mm	BM4010	F	
		GIO 3.3 11	3.3 mm	11 mm	BM4011	F	
		GIO 3.3 13	3.3 mm	13 mm	BM4012	F	
		GIO 3.3 15	3.3 mm	15 mm	BM4013	F	
		GIO 4.1 7	4.1 mm	7 mm	BM4014	F	
		GIO 4.1 9	4.1 mm	9 mm	BM4015	F	
		GIO 4.1 11	4.1 mm	11 mm	BM4016	F	
		GIO 4.1 13	4.1 mm	13 mm	BM4017	F	
		GIO 4.1 15	4.1 mm	15 mm	BM4018	F	
		GIO 4.1 17	4.1 mm	17 mm	BM4019	F	
	a) Head Ø	4.8 mm	GIO 4.8 7	4.8 mm	7 mm	BM4020	F
	b) Aesthetic head	1.8 mm	GIO 4.8 9	4.8 mm	9 mm	BM4021	F
	c) Enossal length	7 - 17 mm	GIO 4.8 11	4.8 mm	11 mm	BM4022	F
	d) Enossal Ø	3.3 / 4.1 / 4.8 mm	GIO 4.8 13	4.8 mm	13 mm	BM4023	F



Delivery incl. surgical screw CST GTI

NOTE IN CONTRAINDICATIONS implants made from c.p. titanium are in general less force resistant compared to those made from titanium alloy. implants having a nominal diameter of less than 3.8 mm are not indicated for single tooth replacement, independently if they are made from c.p. titanium or from titanium alloy. implants Typ GCO 3.7 and GIO 3.3 should be used as supporting implants.

GCO 3.7 mm and GIO 3.3 mm implants may not be used in the molar region and not as single tooth implants with off axis load. GCO 3.7 mm and GIO 3.3 mm implants are used as supporting implants, for example to increase the number of abutments in immediate load situation. Never use GCO 3.7 mm and GIO 3.3 mm implants when cantilevers are involved. GCO 3.7 and GIO 3.3 are not for single tooth replacement.


Note the following contra-indications:

GCO 3.7, GIO 3.3: Do not use in areas where off-axis load is present. Do not use in areas with strong chewing forces. Do not use in flexion areas of the jaws.

GCO IMPLANTS



Precision conical implant with octagon socket and 8° inner cone. Octagonal system. Connection type: Inner cone + internal octagon / head-Ø: 4.8 mm. Strengthened titanium alloy (Ti6Al4V) according to ASTM F 136/ISO 5832-3 or Ti6 Al7 Nb. Compression screw thread and internal cone 8° for stable endosseous anchorage





- a) Head Ø 4.8 mm
- b) Aesthetic head 1.8 mm
- c) Enossal length 7 - 15 mm
- d) Enossal Ø 3.7 / 4.1 / 4.8 mm

Description	Enossal Ø	Enossal length	REF	Price cat.
GCO 3.7 9	3.7 mm	9 mm	BM6510	F
GCO 3.7 11	3.7 mm	11 mm	BM6511	F
GCO 3.7 13	3.7 mm	13 mm	BM6512	F
GCO 3.7 15	3.7 mm	15 mm	BM6513	F
GCO 4.1 7	4.1 mm	7 mm	BM6514	F
GCO 4.1 9	4.1 mm	9 mm	BM6515	F
GCO 4.1 11	4.1 mm	11 mm	BM6516	F
GCO 4.1 13	4.1 mm	13 mm	BM6517	F
GCO 4.1 15	4.1 mm	15 mm	BM6518	F
GCO 4.8 9	4.8 mm	9 mm	BM6524	F
GCO 4.8 11	4.8 mm	11 mm	BM6525	F
GCO 4.8 13	4.8 mm	13 mm	BM6526	F
GCO 4.8 15	4.8 mm	15 mm	BM6527	F




Delivery incl. surgical screw CST STI

SURGICAL ACCESSORIES FOR GIO, GCO

Description	Code	REF	Price cat.
 Cover screw. Tighten with HT 1.25 For 2 mm gingival height, cylindrical	CST GTI	BM6554	A
	HS2 GTI	BM6549	B
 Gingivaformer Tighten with HT 1.25 For 4 mm gingival height, cylindrical	HS4 GTI	BM6550	B
	HS5 GTI	BM6551	B







TEMPBASE® ABUTMENTS FOR PROVISIONAL RESTORATIONS


Temporary basis for GIO® and GCO®. With anti-rotation, PEEK, individually grindable. For cemented crowns and bridges. Tighten with **HT 1.25**.

	TempBase® for GIO, GCO, 2 mm gingival height, incl. screw	TPB 2 GIO	BM6547	D
	TempBase® for GIO, GCO, 4 mm gingival height, incl. screw	TPB 4 GIO	BM6548	D

ABUTMENTS FOR CEMENTED PROSTHETICS FOR GIO, GCO

With octagon. Anti-rotation, screw-in abutment with pre-assembled screw for cemented crowns or bridges.. Tighten with **TT 1.25**.

	Description	Code	REF	Price cat.
	Straight abutment, 6 mm height above implant	TLA GCO	BM4042	D
	15° angled. The angulation runs axial above the edge of the octagon	TLA15 GCO	BM4040	E
	15° angled. The angulation runs axial above the area of the octagon Metalic-yellow	TLA15F GCO	BM4050	E
	20° angled. The angulation runs axial above the edge of the octagon	TLA20 GCO	BM3177	E
	25° angled. The angulation runs axial above the edge of the octagon	TLA25 GCO	BM4041	E
	25° angled. Die angulation runs axial above the area of the octagon Metalic-yellow	TLA25F GCO	BM4069	E



Description	Impression post With Pick-Up screw	Lab analogue Inner octagon	Synthetic base Castable 7.4 mm high For TLA GCO	Synthetic base Castable 4 mm high For TLA15 GCO	Synthetic base Castable 4 mm high For TLA20 GCO	Synthetic base Castable 4 mm high For TLA25 GCO
Code	OLT GCO	IA GCO	PS	PS 15	PS 20	PS 25
REF	BM6571	BM6572	BM6575	BM6576	BM6577	BM6578
Price cat.	B	B	A	A	A	A

Description	Unit	Code	REF	Price cat.
Replacement screw for OLT GCO, long	1 piece	SF 767	BM6661	B

ALTERNATIVE Screw SF OLT, short, 14.5 mm for OLT GCO 1 piece SF OLT BM6579 B

With cone (without anti-rotation). One piece abument with solid thread.



Description

One-piece abutment with solid thread, for cemented crowns or bridges, height above implant 6.0 mm, 8° incline. Can be shortened or trimmed. Flat on one side. Direct impression or transfer to lab analogue IA STI. Tighten with **HT 1.25**.

Code

TCA GTI

REF

BM6594

Price cat.

D



Description

Impression post
9 mm high

Impression post
13 mm high

Lab analogue
With conus

Plastic base,
Castable, 7.4 mm

Code

TS GTI

TSX GTI

IA GTI

PS

REF

BM6599

BM6600

BM6601

BM6575

Price cat.

B

B

B

A

SOLID ABUTMENTS



Description

Flat on one side for increased rotation stability, 6° Incline
Tighten with **IT TCA**

Height above implant

4 mm

5.5 mm

7 mm

Colour

yellow

grey

blue

Code

SA4 GTI

SA5 GTI

SA7 GTI

REF

BM6617

BM6618

BM4044


Price cat.

C

C

C

Tightening torque 30 Ncm. For the impression taking on the abutment use:



Description

Colour coded transfer post for SA-abutments

Use with **REF 460230**

Abutment analogue

Use with **TZ SA4, 5, 7**

Castable abutment
10 mm high, reducible

Anti-rotation (PA SA)
Internally round (PA SR)

Gingiva retractor,
Pack of 4

Use instead of retraction cord

Code

TZ SA4, 5, 7

AA SA4, 5, 7

PA SA/ PA SR

REF

BM6620, BM6621, BM6622

BM6623, BM6624, BM6625

BM6626, BM6627

BM1319

Price cat.

A


B

A

B

Endosseous dental implant systems **GIO, GCO and GIW**

With cone (without anti-rotation). Single-piece abutments with fixed thread and vertical cement escape grooves.



	Description	Height above implant	Colour	Code	REF	Price cat.
	6° incline Tighten with IT TCA	4 mm	yellow	TCA4 GTI	BM6628	C
		5.5 mm	grey	TCA5 GTI	BM6629	C
		7 mm	blue	TCA7 GTI	BM6630	C

Tightening torque 30 Ncm. The impression taking is carried out on the abutment directly.



Description	Abutment analogue Use with TCA 4	Abutment analogue Use with TCA 5	Abutment analogue Use with TCA 7	Castable abutment 10 mm high Reducible, white
Code	AA TCA 4	AA TCA 5	AA TCA 7	PA TCA
REF	BM1101	BM1102	BM1321	BM4030
Price cat.	B	B	B	A

ABUTMENTS FOR SCREW-IN PROSTHETICS

	Description	Code	REF	Price cat.
	With outer octagon and octagon Abutment with double threaded screw and octagon for crowns, bridges and bases. Anti-rotation anchorage of the prosthetics and to the implant. Height above implant 1.5 mm. Tighten with HT 1.25	OSA GCO	BM6580	E
	With outer octagon and conus Singel piece , octagonal screw-in abutment for anti-rotational anchorage of the prosthetics. For crowns, bridges and bar connectors. Impression taking must be taken on the abutment. Height above implant 1.5 mm. Tighten with HT 1.77 / HTX 1.77 . Tightening torque 30-35 Ncm	OSA GTI	BM6650	D



Description	Transfer post For OSA GCO/STI	Screw For anchorage of TST GTI	Lab analogue For OSA GCO/STI	Castable abutment 10 mm high, pack of 5	Screw for fixation
Code	TST GTI	SF 365	OA GTI	Internally edged PAOA	Internally round PAOR
REF	BM6581	BM6582	BM6584	BM6586	BM6587
Price cat.	B	B	B	B	B
					SF 350 or SF 365 BM6588 or BM6582

PROTECTION CAP FOR OSA STO AND OSA STI**Description**

Protection cap made of plastic, white
Screwable and reducible
For OSA GCO or OSA GTI, 3 mm high

Unit

Pack of 5

Code

CGCO

REF

BM6636

Price cat.

B

BAR COMPONENTS**Description**

Bar coping made of CoCrMo for all S-implants
Height above implant 5.5 mm
Screw directly onto implant with SF 350 or SF 365

REF

SK C GTI

Price cat.

C



Bar sleeve made of POM, castable, large
3.0 mm high, 2.2 mm wide, 55 mm long

REF

PA SP1

Price cat.

C



Titanium bar matrix for PA SP1, 55 mm long

REF

SP6

Price cat.

C

CONICAL CROWNS AND BRIDGES**Description**

Bur cylinder for bridges (internally round), height above implant 8 mm, 6 mm Ø
Screw onto OSA GTI or OSA GCO with screw SF 350

REF

KKR

Price cat.

B



Bur cylinder for crowns (internal octagon), height above implant 5 mm, 7 mm Ø
Screw onto OSA GTI or OSA GCO with screw SF 350

REF

FZ 5

Price cat.

B



Bur cylinder for crowns (internal octagon), height above implant 7 mm, 7 mm Ø
Screw onto OSA GTI or OSA GCO with screw SF 350

REF

FZ 7

Price cat.

B

TITANIUM BASE FOR CAD CAM

Titanium base for GIO and GCO. Material Ti6Al4V, anti-rotation. Tighten with **TT 1.25**.

**Description**

Height 3.9 mm above implant shoulder
Incl. screw SF N62

Code

MB4 GIO

REF

BM7759

Price cat.

D

RETENTIVE BALL ABUTMENTS



Description

Ball abutment fits all nylon caps
Head Ø 2.5 mm, height above implant 3 mm
Tighten with **HT 1.25**

Ball abutment fits zu nylon caps
Head Ø 2.5 mm, height above implant 3.6 mm
Tighten with **Tool E** (alternatively with **HAS 2**)

Code **REF** **Price cat.**

TB GTI **BM6653** **D**

TB3 GTI **BM6655** **D**

Impression taking with transfer post

Direct impression taking



Description

Lab analogue

Transfer post
Tighten with **HT 1.25**

Ball analog

Code

IA GTI

TS GTI

IAB

REF

BM6601

BM6599

BM7763

Price cat.

B

B

C

ACCESSORIES FOR BALLHEAD ABUTMENTS



Description

Nylon cap transparent, **Pull-off force** ca. 1200g
(EXTERNAL PRODUCT)

Unit

Pack of 2

Code

NC

REF

BM1325

Price cat.

A1



Nylon cap pink, **Pull-off force** ca. 800g
(EXTERNAL PRODUCT)

Pack of 2

NC 1

BM1326

A1



Nylon cap yellow, **Pull-off force** ca. 500g
(EXTERNAL PRODUCT)

Pack of 2

NC 2

BM1327

A1



Green, strong Nylon caps R-NC
With increased friction strength
Only with reduced diameter ball
≤ 2.3 mm

Pack of 2

R-NC

BM3090

A1



Pink, medium
(EXTERNAL PRODUCT)

Pack of 2

R-NC 1

BM3089

A1



Orange, soft

Pack of 2

R-NC 2

BM3088

A1



Metal sleeve for all nylon caps
(EXTERNAL PRODUCT)

H

BM1329

B

RETENTIVE BALL ABUTMENTS



Description

Ball abutment fits gold cap GC
Head Ø 2.3 mm, height above implant 3.4 mm
Tighten with **Tool E** (alternatively with **HAS 2**)

Code

TB2 GTI

REF

BM6654

Price cat.

D

Impression taking with transfer post

Direct impression taking



Description

Lab analogue with conus

Transfer post
Tighten with **HT 1.25**

Ball analog

Gold cap

Code

IA GTI

TS GTI

IAB 2

GC

REF

BM6601

BM6599

BM7771

BM1106

Price cat.

B

B

C

upon request

LOCALICER®

Abutment for removable prosthetics. We recommend a minimum of six implants per jaw and the use of a single denture as splint when using LOC abutments. Tighten with **HT 1.77**.



Description

LOC S 2.5

LOC S 4

Height

2.5 mm

4.0 mm

REF

BM6592

BM6593

Price cat.

D

D

ACCESSORIES FOR LOCALICER®



Description

Analogue + impression set

Set with 5 caps + 1 housing (**EXTERNAL PRODUCT**)

Pull-off force

Yellow 600 g, Pink 1.200 g, Transparent 1.800 g, Violet 2.700 g
Black has no retention and is designed for temporary solutions for up to one month

Code

AA LOC

NCS

REF

BM3142

BM3143

Price cat.

C

D


GIW IMPLANTS



OCTAGONAL SYSTEM

Implant with internal octagon and 8° internal cone. For anti-rotation superstructures. Highly fracture resistant titanium alloy (Ti6Al4V) according to ASTM F 136 / ISO 5832-3.

These implants with a head diameter of 6.5 mm can be used, if adequate bone is available, for posterior single crowns with the width of one premolar in dentate and partially/fully edentulous ridges as for superstructures retained by bars and ball attachments.



Description	Enossal Ø	Enossal length	REF	Price cat.
GIW 4.8 9	4.8 mm	9 mm	BM4049	F
GIW 4.8 11	4.8 mm	11 mm	BM4045	F
GIW 4.8 13	4.8 mm	13 mm	BM4046	F

- a) Head Ø 6.5 mm
- b) Aesthetic head 1.8 mm
- c) Enossal length 9 - 13 mm
- d) Enossal Ø 4.8 mm


Delivery incl. surgical screw CST GIW



SURGICAL ACCESSORIES FOR GIW

Description	Code	REF	Price cat.
 Cover screw. Tighten with HT 1.25	CST GIW	BM4051	B
 Gingivaformer Tighten with HT 1.25 For 2 mm gingival height Cylindrical	HS2 GIW	BM4052	B
 Gingivaformer Tighten with HT 1.25 For 4 mm gingival height Cylindrical	HS4 GIW	BM4053	B

TITANIUM BASE FOR CAD CAM

Description	Code	REF	Price cat.
 Titanium base for GIW, Material Ti6Al4V, anti-rotation. Height 3.9 mm above implant shoulder. Incl. screw SF N62 Tighten with TT 1.25	MB4 GIW	BM7762	D

ABUTMENTS FOR CEMENTED PROSTHETICS FOR GIW

Anti-rotation screw-in abutment incl. screw, for cemented crowns and bridges. Tighten with TT 1.25

**Description**

Straight abutment

Code

TLA GIW

REF

BM4054

Price cat.

D

15° angled abutment

The 15° line angle runs axially across the edge of the octagon

TLA15 GIW

BM4055

E

**Description**Impression post with
Pick-Up screw
Tighten with HT 1.25

Lab analogue

Plastic base
Castable
For TLA GIWPlastic base
Castable
For TLA15 GIW**Code**

OLT GIW

IA GIW

PA GIW

PA15 GIW

REF

BM4056

BM4057

BM4059

BM4060

Price cat.

B

B

A

A

Description

Replacement screw for TLAs of all S-implants

Code

SF TLA

REF

BM6598

Price cat.

B

Replacement screw for impression post OLT GIW

SF 767

BM6661

B

ALTERNATIVE

Screw SF OLT, short, 14.5 mm for OLT GIW

SF OLT

BM6579

B

SOLID ABUTMENTS**Description**Flattened on one side for increased
rotational stability
6° angled
Tighten with HT 1.77**Height above implant**

4 mm

6 mm

Code

SA4 GIW

SA6 GIW

REF

BM4061

BM4062

Price cat.

C

C

**Description**Transfer post for
SA6 abutmentAbutment analogue
Brown anodized
For TZ SA6 GIWPlastic base
Castable**Code**

TZ SA6 GIW

AA SA6 GIW

PA GIW

REF

BM4063

BM4064

BM4059



Price cat.


A

B

A








RETENTIVE BALL ATTACHMENTS

	Description	Code	REF	Price cat.
	Ball abutment, fits all nylon caps Head Ø 2.5 mm Height above implantat 3 mm Tighten with HT 1.25	TB GIW	BM4068	D
	Ball abutment, fits all nylon caps Head Ø 2.5 mm Height above implantat 3.45 mm Tighten with Tool E (alternatively with HAS 2)	TB3 GIW	BM4066	D







Description	Lab analogue with conus	Transfer post Tighten with HT 1.25	Ball analog
Code	IA GIW	TS GIW	IAB
REF	BM4057	BM4065	BM7763
Price cat.	B	B	C




ACCESSORIES FOR BALLHEAD ABUTMENTS

	Description	Unit	Code	REF	Price cat.
	Nylon cap transparent, Pull-off force ca. 1200g (EXTERNAL PRODUCT)	Pack of 2	NC	BM1325	A1
	Nylon cap pink, Pull-off force ca. 800g (EXTERNAL PRODUCT)	Pack of 2	NC 1	BM1326	A1
	Nylon cap yellow, Pull-off force ca. 500g (EXTERNAL PRODUCT)	Pack of 2	NC 2	BM1327	A1
	Green, strong Nylon caps R-NC With increased friction strength Only with reduced diameter ball ≤ 2.3 mm	Pack of 2	R-NC	BM3090	A1
	Pink, medium (EXTERNAL PRODUCT)	Pack of 2	R-NC 1	BM3089	A1
	Orange, soft	Pack of 2	R-NC 2	BM3088	A1
	Metal sleeve for all nylon caps (EXTERNAL PRODUCT)		H	BM1329	B

DRILL SEQUENCES CYLINDRICAL IMPLANTS

Implants	Enossal Ø	recommended drill sequences
 GIO	3.3 mm	DS 2 BM1359 → DS 2.2 BM1405 → DS 2.8 BM1404 → DRA 3.3 CS BM1378
		DSL 2 BM6555 → DSL 2.2 BM6556 → DSL 2.8 BM6557 → DRA 3.3 CS BM1378
 GIO	4.1 mm	DS 2 BM1359 → DS 2.2 BM1405 → DS 2.8 BM1404 → DS 3.5 BM1375 → DRA 4.1 CS BM1379
		DSL 2 BM6555 → DSL 2.2 BM6556 → DSL 2.8 BM6557 → (DLS 3.5) (BM6559) → DRA 4.1 CS BM1379
 GIO	4.8 mm	DS 2 BM1359 → DS 2.2 BM1405 → DS 2.8 BM1404 → DS 3.5 BM1375 → DS 4.0 BM1376 → (DS 4.2) (BM1377) → DRA 4.1 CS BM1379
		DSL 2 BM6555 → DSL 2.2 BM6556 → DSL 2.8 BM6557 → DLS 3.5 BM6559 → DSL 4.0 BM6560 → (DSL 4.2) (BM6561) → DRA 4.1 CS BM1379
 GIW	4.8 mm	DS 2 BM1359 → DS 2.2 BM1405 → DS 2.8 BM1404 → DS 3.5 BM1375 → DS 4.0 BM1376 → (DSL 4.2) (BM6561) → DRA 4.8 CSW BM6568
		DSL 2 BM6555 → DSL 2.2 BM6556 → DSL 2.8 BM6557 → DSL 3.5 BM6559 → DSL 4.0 BM6560 → (DSL 4.2) (BM6561) → DRA 4.8 CSW BM6568

DRILL SEQUENCES CONICAL IMPLANTS


Implants	Enossal Ø	Length in mm	recommended drill sequences
 GCO	3.7 mm	9 / 11 / 13	DS 2 BM1359 → DS 2.2 BM1405 → DRK 1 BM6562 → DRA 3.3 CS BM1378
			DSL 2 BM6555 → DSL 2.2 BM6556 → DRK 1 BM6562 → DRA 3.3 CS BM1378
	15		DS 2 BM1359 → DS 2.2 BM1405 → DRK 2 BM6563 → DRA 3.3 CS BM1378
			DSL 2 BM6555 → DSL 2.2 BM6556 → DRK 2 BM6563 → DRA 3.3 CS BM1378
 GCO	4.1 mm	7 / 9 / 11 / 13	DS 2 BM1359 → DS 2.2 BM1405 → DRK 3 BM6564 → DRA 4.1 CS BM1379
			DSL 2 BM6555 → DSL 2.2 BM6556 → DRK 3 BM6564 → DRA 4.1 CS BM1379
	15		DS 2 BM1359 → DS 2.2 BM1405 → DRK 2 BM6563 → DRA 4.1 CS BM1379
			DSL 2 BM6555 → DSL 2.2 BM6556 → DRK 2 BM6563 → DRA 4.1 CS BM1379
 GCO	4.8 mm	7 / 9 / 11	DS 2 BM1359 → DS 2.2 BM1405 → DRK 4 BM6565 → DRA 4.1 CS BM1379
			DSL 2 BM6555 → DSL 2.2 BM6556 → DRK 4 BM6565 → DRA 4.1 CS BM1379
	13 / 15		DS 2 BM1359 → DS 2.2 BM1405 → DRK 5 BM6566 → DRA 4.1 CS BM1379
			DSL 2 BM6555 → DSL 2.2 BM6556 → DRK 5 BM6566 → DRA 4.1 CS BM1379

PLEASE NOTE that GCO is by design a compression screw and a drill sequence can not be recommended for all bone qualities. In low density bone a smaller drill is recommended.

GUIDE JACKET




Description	Unit	Material	REF	Price cat.
BFH 2.0 guide jacket 2.0 mmd	Pack of 5	Ti6Al4V	BM7100	B
BFH 2.5 guide jacket 2.5 mmd	Pack of 5	Ti6Al4V	BM7101	B
BFH 3.0 guide jacket 3.0 mmd	Pack of 5	Ti6Al4V	BM7102	B
BFH 3.2 guide jacket 3.2 mmd	Pack of 5	Ti6Al4V	BM7103	B
BFH 3.5 guide jacket 3.5 mmd	Pack of 5	Ti6Al4V	BM7104	B



a) Length	5 mm
b) Height of step	0.7 mm
c) Max. Ø top	3.7 - 5 mm
d) Nominal Ø	3 - 4.4 mm
e) Ø of drilling in the drill template	2.05 - 3.55 mm

PATHFINDER DRILLS

	Description	Type	Working length	REF	Price cat.
	BCD 1	non sterile	15 mm	BM2100	C
	BCDX 1	non sterile	15 mm	BM2103	C

TWIST DRILL (EXTERNALLY IRRIGATED) FOR GIO, GIW







Surgical steel, laser marked. Non sterile. For lifetime extension and consistent cutting performance and easy removal of contaminated areas. Maximum speed with external cooling 700 rpm. Drills can be used without cooling at 50 rpm on green contra-angle handgrip.

These drills are made of premium material using state-of-the-art precision milling machines. The benefit for you and your patients: Extremely good and durable cutting performance due to sophisticated geometry of the blades.

	Description	Length	Working length	Drill Ø	REF	Price cat.
	DS 2	32.5 mm	17 mm	2 mm	BM1359	D
	DS 2.2	36.5 mm	18 mm	2.2 mm	BM1405	D
	DS 2.8	36.5 mm	18 mm	2.8 mm	BM1404	D
	DS 3.2	36.5 mm	18 mm	3.2 mm	BM1403	D
	DS 3.5	36.5 mm	18 mm	3.5 mm	BM1375	D
	DS 4.0	36.5 mm	18 mm	4.0 mm	BM1376	D
	DS 4.2	36.5 mm	18 mm	4.2 mm	BM1377	D
	DSL 2	45.5 mm	27 mm	2 mm	BM6555	D
	DSL 2.2	45.5 mm	27 mm	2.2 mm	BM6556	D
	DSL 2.8	45.5 mm	27 mm	2.8 mm	BM6557	D
	DSL 3.2	45.5 mm	27 mm	3.2 mm	BM6558	D
	DSL 3.5	45.5 mm	27 mm	3.5 mm	BM6559	D
	DSL 4.0	45.5 mm	27 mm	4.0 mm	BM6560	D
	DSL 4.2	45.5 mm	27 mm	4.2 mm	BM6561	D

-55%
heat**HEATLESS® DRILLS DRK FOR GCO IMPLANTS**




Form drills made from surgical steel (SS), length and color coded, resterilizable. Note the information in the **Instructions for Use** for hygiene and preparation requirements with medical products intended for multiple use in implantology. Pilot drilling using DS2, DSL2 or Pathfinder (see previous page).

	Description	System	Colour	Length	Implantat 0	REF	Price cat.
	DRK 1	GCO	yellow	9, 11, 13 mm	3.7 mm	BM6562	D
	DRK 2	GCO	black	15, 17, 19 mm	3.7 / 4.1 mm	BM6563	D
	DRK 3	GCO	red	7, 9, 11, 13 mm	4.1 mm	BM6564	D
	DRK 4	GCO	blue	7, 9, 11 mm	4.8 mm	BM6565	D
	DRK 5	GCO	green	13, 15 mm	4.8 mm	BM6566	D
	DRK 6	GCO	metallic	21, 23, 25 mm	3.7 / 4.1 mm	BM6567	D






IT HAS BEEN SCIENTIFICALLY PROVEN

Heatless® drills by Dr. Ihde Dental generate 55 % less heat than traditional bone drills from other manufacturers. This makes it possible to use higher rotational speeds: between 3,000 and 5,000 rpm are recommended with good external cooling and intermittent drill technique.

COUNTERSINK HEADS

	Description	System	Type	Implantat 0	REF	Price cat.
	DRA 3.3 CS	GIO	non sterile	3.3 mm	BM1378	D
	DRA 4.1 CS	GIO, GCO	non sterile	4.1 mm	BM1379	D
	DRA 4.8 CSW	GIW	non sterile	4.8 mm	BM6568	D

TAP (MADE OF CRONIDUR) FOR GIO, GIW

	Description	Finish	Length	Implantat 0	REF	Price cat.
	TAP1 3.3	9-17 mm	34 mm	3.3 mm	BM1381	D
	TAP1 4.1	9-17 mm	34 mm	4.1 mm	BM1382	D
	TAP1 4.8	9-17 mm	30 mm	4.8 mm	BM1383	D
	TAP2 3.3	9-17 mm	23 mm	3.3 mm	BM6569	D
	TAP2 4.1	9-17 mm	23 mm	4.1 mm	BM6570	D


DRILLSTOP TRAY

Autoclaveable up to 134° C, not suitable for dry heat sterilizers.











Description	REF	Price €
Drillstop A	BM1500	
Drillstop C	BM1502	
Drillstop F	BM1505	
Drillstop H	BM1507	
Drillstop K	BM1510	
Drillstop L	BM1511	
Drillstop M	BM1512	
Drillstop N	BM1513	
Drillstop O	BM1514	
Drill DS 2.8	BM1404	
Drill DS 3.2	BM1403	
Drill DS 3.5	BM1375	
Drill DS 4.0	BM1376	
Drill DS 4.2	BM1377	
Drill DSL 2.8	BM6557	
Drill DSL 3.2	BM6558	
Drill DSL 3.5	BM6559	
Drill DSL 4.0	BM6560	
Drill DSL 4.2	BM6561	
Tray with content	BM6677	779.00

INSERTION TOOLS FOR GIO, GCO, GIW

	Description	Type	Length	Index	For	REF	Price cat.
	IT1 GCO	long	20 mm	HT 1.25		BM1384	K
	IT2 GCO	short	12 mm	HT 1.25		BM1385	K
	IT3 GCO	medium	16 mm	HT 1.25		BM1386	K
	ITW GIO	contra-angle	23 mm			BM1389	F
	ITWH GIO	contra-angle / hex	23 mm			BM1391	F
	This instrument is with an additional hex for force transmission to the contra-angle handpiece equipped. Fits to W&H contra-angle handpiece only.						
	ITV	short	11 mm		GIO, GCO	BM1390	C
	IT ITV	Ratchet adapter			Adapter for ITV	BM1365	D
	Adapter	short / contra-angle	22 mm		ITV 500850	BM1366	D
	Adapter	long / contra-angle	32 mm		ITV 500850	BM1367	D
	Adapter	medium / contra-angle	27 mm		ITV 500850	BM1368	D
	HAS	square	85 mm		ITV	BM1373	H
	IT TCA	long	20 mm	direct	SA GTI, TCA GTI	BM1322	D
	Tool E	long	20 mm	direct	TB2 GTI, TB3 GTI, TB3 GIW	BM3336	D

INSTRUMENTS FOR SCREWS

	Description	Type	Length	Ø	Code	REF	Price cat.
	Hex instrument	long	21 mm	1.25 mm	HT 1.25	BM3022	C
	Hex instrument	long, for contra-angle	26 mm	1.25 mm	HT 1.25 M	BM3047	C
	Hex instrument	short	14 mm	1.25 mm	HTS 1.25	BM3023	C
	Hex instrument	extralong	45 mm	1.25 mm	HTX 1.25	BM7764	C
	Hex instrument	long	19 mm	1.77 mm	HT 1.77	BM3024	C
	Hex instrument	extralong	45 mm	1.77 mm	HTX 1.77	BM1070	C
	Hex instrument	long, for contra-angle	28 mm	1.77 mm	HT 1.77 M	BM3048	C
	Torx-Instrument	for GIO, GCO and GIW abutments	21 mm	1.25 mm	TT 1.25	BM3027	C

TOOLS

	Description	Type	Code	REF	Price cat.
	Guide jacket	For pilot drill, Titanium 10 mm, 2.2 mm Ø. Pack of 5	BFH	BM6537	B
	X-ray measuring sphere	Surgical steel, Ø 5 mm Pack of 5	RM	BM3006	A
	Drill extension contra-angle	Extends by 19 mm	DX2	BM1349	D
	Punch	For contra-angle, 4.9 mm Ø	PUW1	BM3002	C
	Standardized probe	Scale 1 mm for radiographical measuring for cylindrical implant types, 22 mm	PDG	BM1350	A
	HAS	Flat spanner for ITV insertion tool, square-end, 85 mm	HAS	BM1373	H
	HAS 2	Flat spanner for ball attachments TB2 GT1, TB3 GT1 and TB3 GIW, 85 mm	HAS 2	BM1392	H
	Ratchet	For all Hex instruments and insertion tools, without torque	RAT 2	BM1352	K
	Torque wrench	Torque wrench 10 - 70 Ncm for all insertion tools, hex- and torx instruments	TW2	BM1356	S
	Universal adapter	For all contra-angle instruments Use with ratchet RAT 2. Max. 30 Ncm	UAW	BM3026	E

STARTER TRAY

Autoclaveable up to 134° C, not suitable for dry heat sterilizers. This surgical kit contains all drills and tools for first works with the S-System.
Material: autoclavable plastic.



Description	Code	REF	Price €
Insertion tool for ratchet	IT3 GCO	BM1386	
Insertion tool for contra angle	ITWH GIO	BM1391	
Hex-instrument	HT 1.25	BM3022	
Torx instrument	TT 1.25	BM3027	
Pilot drill	BCD 1	BM2100	
Form drill	DS 2.2	BM1405	
Form drill	DS 2.8	BM1404	
Form drill	DS 3.5	BM1375	
Form drill	DS 4.2	BM1377	
Countersink head	DRA 4.1 CS	BM1379	
Standardized probe	PDG	BM1350	
Standardized probe	PDG	BM1350	
Torque wrench	TW2	BM1356	
Starter tray with content		BM6501	upon request
Starter tray w/o content		SBM4009	upon request

ONEWAY BIOMEED

We are certified DIN EN ISO 13485, and annex II of EEC Directive 93/42 EWG (2007).
Due to technical reasons the product dimensions shown in this brochure might deviate from reality.
GIO, GCO and GIW are registered trademarks.

In case that implants would be reprocessed (cleaned, resterilized) infections could occur,
because no validated procedures for reprocessing are available.

(The products of this catalogue are CE marked (class I) and CE 1936 marked (class IIa and IIb) according to 93/42/EC Directive).

Commercial products that are not monitored by our notified body are declared as third-party products.

Compilation and explanation of symbols on the packaging:



Batch No.



Sterilized by radiation



Non-sterile



Intended for use by dentists or surgeons only



Single use product



Instruction for use



Expiry date



Store in a dry place



Store tightly keep closed



Do not use if packing is damaged



Do not resterilize



Manufacturer



Production date



Catalogue number



Secure rotation protection by
precision inner octagon and
internal cone 8°

Reliable insertion - ease of use

Universal application for
permanent and removable
prosthodontics

ONEWAY
BIOMEED



Dr. Ihde Dental AG
Dorfplatz 11
CH - 8737 Gommiswald / SG
Tel +41 (0)55 293 23 23
contact@implant.com
www.implant.com

Distributed by **onewaybiomed GmbH**
Dorfplatz 11
CH - 8737 Gommiswald / SG
Phone +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.onewaybiomed.com



Dr. Ihde Dental GmbH
Erfurter Str. 19
D - 85386 Eching / München
Tel +49 (0)89 319 761 0
info@ihde-dental.de