

ONEWAY
BIOMED

IMMEDIATE LOADING single part implants

DENTAL IMPLANT SYSTEM GBL[®]

SYMBOLS FOR IMPLANT PROPERTIES AND PROSTHETIC SOLUTIONS

Small Abutment Head



Multi-Unit Abutment



Apical Wide Thread



Large Abutment Head



Micro Thread



For Screw-in Prosthetics




Ball Abutment



Bendable



For Cemented Prosthetics



Localicer® Abutment



ONEWAY
BIOMED

THE ADVANTAGES

OF THE ENDOSSEOUS DENTAL IMPLANT SYSTEM GBL®

The surface of GBL® implants provide a specially lasered surface with exactly defined properties. For anti-rotation an internal square connects with press-fit to the abutment. The cone in combination with the internal stare provides stability and 100% tightness. GBL® implants are universally applicable for fixed and removable prosthetics.

The prescribed or recommended tightening torques for implants, abutments and screws can be found on our website:

www.implant.com/en/downloads




Safely
antirotational
thanks to its internal
precision square

Cone technology
for a tight seal

Universally
suitable for fixed
and removable
prosthodontics

Made of highly
resistant
titanium alloy

Smart
instrument tray




PREPARATORY STEPS WHEN USING A DRILLING TEMPLATE

1. Have your laboratory produce a drilling template with the appropriate drill holes for the marker bore. To be on the safe side, the laboratory might insert guide sleeves (**REF BFH**) into the drill holes to ensure that the drilling angle is correct. Use a 2 mm \varnothing drill for pilot drilling.
2. For subsequent drilling sequences, drill stops can be used that are slid over the drill according to the appropriate depth of the drill hole and screwed in place. Consider the thickness of the mucosa and the height of the template as appropriate.

Thanks to the extremely high cutting efficiency of our drills, no ascending drilling sequences will usually be required.


Recommended RPM: 2000-5000

Apply sufficient cooling and allow the cooling to reach the working blades of the drills. Drill stop taking from \varnothing 2.8.



SURGERY

1. Recommended drill sequence



Owing to the high quality and geometry of the blades of our drills, the final preparation may be performed immediately after the pilot drilling.


2. Implant packaging



Original packaging




Open the sealed cover at the lid.
Remove the label and place it into
the patients record.



The open pack contains the
implant in a sterile tube
(primary packaging).

3. Remove the implant from its packaging

1. Open the Lid
2. The implant is attached to the cap and can be removed by breaking it off at the pre-determined breaking point
3. Remove the implant, making sure not to touch the inner wall of the tube





4. Handling

4.1 Connect

Attach the placement aid to the implant, holding the cap to which the implant is attached with the other hand.

4.2 Mounting the adapter ITV WST / contra-angle


Place the ITV Wst (angled handpiece) or IT ITV (ratchet) adapter on the ITV GBL placement aid. Mount the placement aid. Hold the cap firmly in one hand and break off the implant at the pre-determined breaking line.



4.3 Alternative to 4.2:

Place the IT ITV (ratchet) adapter on the ITV GBL placement aid.

Mount the placement aid.
Hold the cap firmly in one hand and break off the implant at the pre-determined breaking line.




5. Insertion

Use the angled handpiece, ratchet or shank to screw the implant into the implant bed (clockwise).

The enossal aspect of the implant must be submerged in the bone.
Upon **complete** insertion, the implant may be turned back ¼ revolution to reduce the load on the bone.

The system is suitable for deep insertion (below bone level).



6. Remove insertion tool from implant

Loosen the insertion aid from the implant by pulling it off.




7. Result



8. Aftercare

Seal the implant with a matching cover screw.



After healing time:
Remove surgical screw.

Arrow: Top view.
Rotation counter clockwise.




9. Pick-up impressions

9.1 Impression with perforated custom tray

Torx-instrument TT 1.25

Insert impression posts
HLT GBL NC/RC

GBL® implant




9.2 Before taking the impression

Take an impression in an A silicone. You can use the open-tray or the closed-tray technique.

It is necessary to remove the HLT GBL NC/RC impression post from the implant to be able to take out the impression tray.

Impression post HLT GBL NC/RC

GBL® implant



9.3 Taking the impression


Detach the HLT GBL NC/RC from the implant.
HLT GBL NC/RC will stay within the impression.

Use TT 1.25 to loosen screw

Relief window in the
impression tray

HLT GBL NC/RC

GBL® implant




9.4

View of the impression post in the impression
(pick-up technique, bottom view).

Position of the Impression post

Impression material




9.5

Once the impression has been taken, the implant is closed with a healing abutment, while the impression is sent to the laboratory.

TT 1.25

Surgical screw CSTB NC/RC

GBL® implant



10. Closed tray impression taking


10.1 Impression with closed tray

Impression with custom tray.

Securing the impression post
with the thumbscrew

TS GBL NC/RC


GBL® implant



10.2 Before taking the impression

Take an impression in an A silicone.
You can use the open-tray or the closed-tray
technique.


With the closed impression technique, the TS
GBL NC/RC will always remain on the implant
when removing the impression.



10.3 Removing the impression

In the case of closed impressions, the TS GBL
NC/RC impression post will remain on the im-
plant after removing the impression tray.

The impression post will be removed after-
wards.




10.4

Once the impression has been taken, the
implant is closed with an HA NC/RC healing
abutment, while the impression is sent to the
laboratory.

TT 1.25
insert surgical screw
CSTB NC/RC

GBL® implant



11. Procedures in the laboratory


11.1 Pick-up technique

Tighten the IAB against the HLT GBL (NC/RC) impression post.

Use the TT 1.25 to insert the lab analogue

HLT GBL NC/RC

IAB NC or IAB RC



11.2 Closed technique

Secure the IAB NC/RC against the TS GBL (NC or RC) (A)


Reposition the impression post inside the impression (B)

Pour the impression.

Use the thumbscrew to tighten the impression post on the lab analogue.

TS GBL NC/RC

IAB NC or IAB RC




11.3

Pour the impression in dental stone, then remove the impression posts from the lab analogues.

Lab analogue

Fill with gypsum



11.4

The lab analogue will now be embedded in the gypsum in the correct position.

IAB NC/RC



11.5

Positioning of the screw-retained TLA2 15 GBL RC abutment, determining its optimal position and correct angulation.

NOTE The square end must be inserted completely into the analogue.

TT 1.25

Insert screw

TLA2 15 GBL NC/RC
Watch out for the correct
square end position



Arrow: Top view.
Rotation clockwise

IAB NC/RC



11.6

The correct position of the abutment must be ensured during transfer to the mouth.

TLA2 15 GBL NC/RC



11.7

If multiple angled abutments are used, the laboratory will produce a removable resin splint (e.g. from pattern resin) to facilitate positioning within the mouth.


TLA2 15 GBL RC

Pattern Resin




BONE LEVEL PLUS® IMPLANTS

With No-Itis® Laser Surface.
The implant body is made of Ti6Al4V.



		Description	Enossal Ø	Enossal length	REF	Price cat.
		GBL 3.3 8 NC	3.3 mm	8 mm	BM6010	H
		GBL 3.3 10 NC	3.3 mm	10 mm	BM6011	H
		GBL 3.3 12 NC	3.3 mm	12 mm	BM6012	H
		GBL 3.3 14 NC	3.3 mm	14 mm	BM6013	H
		GBL 4.1 8 RC	4.1 mm	8 mm	BM6014	H
		GBL 4.1 10 RC	4.1 mm	10 mm	BM6015	H
		GBL 4.1 12 RC	4.1 mm	12 mm	BM6016	H
		GBL 4.1 14 RC	4.1 mm	14 mm	BM6017	H
a) Enossal Ø	3.3 - 4.8 mm	GBL 4.8 8 RC	4.8 mm	8 mm	BM6018	H
b) Enossal length	8 - 14 mm	GBL 4.8 10 RC	4.8 mm	10 mm	BM6019	H
NC	Narrow Collar	GBL 4.8 12 RC	4.8 mm	12 mm	BM6020	H
RC	Regular Collar	GBL 4.8 14 RC	4.8 mm	14 mm	BM6021	H

Min. Insertion torque 35 Ncm



Square drive

Delivery inclusive insertion tool ITV GBL and surgical screw REF BM6056 or BM6057















- Safely anti-rotational thanks to its internal precision square
- Cone technology for a tight seal
- Universally suitable for fixed and removable prosthodontics
- The cone centers the abutment and provides 100% tightness

SURGICAL SCREWS






Description	Code	REF	Price cat.
Surgical screw for GBL 3.3	CSTB NC	BM6058	B
Surgical screw for GBL 4.1 and 4.8	CSTB RC	BM6059	B




GINGIVAFORMER

		Description	Code	REF	Price cat.
		Gingivaformer conical	GF NC 3.6 2	BM6047	B
		Gingivaformer conical	GF NC 3.6 3.5	BM6048	B
		Gingivaformer conical	GF NC 4.8 3.5	BM6049	B
		Gingivaformer conical	GF RC 4.5 2	BM6050	B
		Gingivaformer conical	GF RC 4.5 4	BM6051	B
		Gingivaformer conical	GF RC 4.5 6	BM6086	B
		Gingivaformer conical	GF RC 6 2	BM6087	B
		Gingivaformer bottle shape	GFB NC 3.3 3.5	BM6052	B
		Gingivaformer bottle shape	GFB NC 3.3 5	BM6053	B
		Gingivaformer bottle shape	GFB RC 4.4 4	BM6054	B
		Gingivaformer bottle shape	GFB RC 4.7 6	BM6055	B

BUR CYLINDER

		Description	Code	REF	Price cat.
		Bur cylinder for GBL 3.3 for telescope crowns	FZB NC	BM6058	D
		Bur cylinder for GBL 4.1 and 4.8 for telescope crowns	FZB RC	BM6059	D
Recommended insertion torque 30 Ncm					

ANALOGUES

		Description	Code	REF	Price cat.
		Implant analogue for GBL 3.3	IA GBL NC	BM6060	B
		Implant analogue for GBL 4.1 and 4.8	IA GBL RC	BM6061	B

STANDARD ABUTMENTS

Description	Code	REF	Price cat.
Abutment for cementing on GBL 3.3, step 1 mm high Height above step 4 mm, incl. screw SFBC NC	CAB 1 NC	BM6062	E
Abutment for cementing on GBL 3.3, step 3 mm high Height above step 4 mm, incl. screw SFBC NC	CAB 3 NC	BM6045	E
Abutment for cementing on GBL 4.1 and 4.8, step 1 mm high Height above step 5.5 mm, incl. screw SFBC RC	CAB 1 RC	BM6063	E
Abutment for cementing on GBL 4.1 and 4.8, step 3 mm high Height above step 5.5 mm, incl. screw SFBC RC	CAB 3 RC	BM6064	E
Recommended insertion torque 20 Ncm			

SCREW-RETAINED ABUTMENTS (REDUCIBLE, GRINDABLE)

Description	Code	REF	Price cat.
Abutment Incl. screw SF B	TAB GBL NC/RC	BM6032	D
Abutment for GBL 3.3, 15° angled Anti-rotational, incl. screw SFB NC	TLA2 15 GBL NC	BM6030	F
Abutment for GBL 4.1 and 4.8, 15° angled Anti-rotational, incl. screw SFB RC	TLA2 15 GBL RC	BM6031	F
Recommended insertion torque 20 Ncm			

ANATOMICAL ABUTMENTS

Description	Code	REF	Price cat.
Anatomical abutment for GBL 3.3 Anti-rotational, incl. screw SFB NC	ANAB NC	BM6067	F
Anatomical abutment for GBL 4.1 and 4.8 Anti-rotational, incl. screw SFB RC	ANAB RC	BM6068	F
Recommended insertion torque 20 Ncm			

TITANIUM BASE FOR CAD CAM**Description**

Titanium base for GBL 3.3, anti-rotation
Incl. screw SFB NC

Titanium base for GBL 4.1 and 4.8, anti-rotation
Incl. screw SFB RC

Code**MB GBL NC****MB GBL RC****REF****BM6076****BM6077****Price cat.****D****D****CASTABLE ABUTMENTS****Description**

Castable abutment for GBL 3.3
Incl. metal base and screw

Castable abutment for GBL 4.1 and 4.8
Incl. metal base and screw

Material

CoCrMo/plastic

CoCrMo/plastic

Code**PLAB2 GBL NC****PLAB2 GBL RC****REF****BM6005****BM6007****Price cat.****G****G****PICK-UP IMPRESSION POST FOR PICK-UP IMPRESSIONS****Description**

Impression post for GBL 3.3

Impression post for GBL 4.1 and 4.8

Code**HLT GBL NC****HLT GBL RC****REF****BM6069****BM6046****Price cat.****C****C****IMPRESSION POST FOR CONVENTIONAL IMPRESSIONS****Description**

Impression post for GBL 3.3

Impression post for GBL 4.1 and 4.8

Impression post long for GBL 3.3

Impression post long for GBL 4.1 and 4.8

Code**TS GBL NC****TS GBL RC****TSL GBL NC****TSL GBL RC****REF****BM6070****BM6071****BM6072****BM6073****Price cat.****C****C****C****C**

ABUTMENTS FOR SCREW-ON PROSTHETIC

	Description	Code	REF	Price cat.
	Gingiva height 0.5 mm	TCT GBL NC 0.5	BM6183	D
	Gingiva height 1.5 mm	TCT GBL NC 1.5	BM6184	D
	Gingiva height 3.5 mm	TCT GBL NC 3.5	BM6185	D
	Gingiva height 0.5 mm	TCT GBL RC 0.5	BM6180	D
	Gingiva height 1.5 mm	TCT GBL RC 1.5	BM6181	D
	Gingiva height 3.5 mm	TCT GBL RC 3.5	BM6182	D
Tighten with HT 1.77				

IMPRESSION TAKING AND LABORATORY ACCESSORIES

In this approach the position of the TCT hex is assigned.

						
	Transfer post	Long screw Tighten with HT 1.25	TCT analogue	Castable abutment 12mm high Internally round Pack of 5	Castable abutment 12mm high Internally edged Pack of 5	Fastening screw Tighten with HT 1.25
Code	TST	SFL	BTT	PSTR (grey)	PSTA	SF
REF	BM3124	BM3125	BM3126	BM3127	BM3128	BM3120
Price cat.	B	B	B	B	B	B

LOCALICER® FOR REMOVABLE PROSTHETIC

If LOC abutments are used in the upper jaw, we recommend to place at least six implants and to splint them through prosthetics in a stable manner. Tighten with HT 1.77.



Description	Height	Code	REF	Price cat.
Localicer® for GBL 3.3	2 mm	LOC GBL NC 2	BM6074	D
Localicer® for GBL 3.3	3 mm	LOC GBL NC 3	BM6075	D
Localicer® for GBL 3.3	4 mm	LOC GBL NC 4	BM6080	D
Localicer® for GBL 4.1 and 4.8	2 mm	LOC GBL RC 2	BM6081	D
Localicer® for GBL 4.1 and 4.8	3 mm	LOC GBL RC 3	BM6082	D
Localicer® for GBL 4.1 and 4.8	4 mm	LOC GBL RC 4	BM6083	D


ACCESSORIES FOR LOCALICER®

Description	Code	REF	Price cat.
Analogue + impression set	AA LOC	BM3142	C
Set with 5 caps + 1 housing (EXTERNAL PRODUCT)	NCS	BM3143	D

Pull-off force
 Yellow 600 g, Pink 1.200 g, Transparent 1.800 g, Violet 2.700 g
 Black has no retention and is designed for temporary solutions for up to one month

MULTI-UNIT ABUTMENTS

Insertion of the angled MU2 abutments with HT 1.25. Insertion of the straight MU2S abutments with HT 1.77














Description	Material	Code	REF	Price cat.
Abutment 17° angled Incl. screw SFB RC	Ti6Al4V	MU2 17 GBL RC	BM6670	L
Abutment 35° angled Incl. screw SFB RC	Ti6Al4V	MU2 35 GBL RC	BM6671	L
Abutment straight Gingiva height 0.5 mm	Ti6Al4V	MU2S 0.5 GBL RC	BM6672	G
Abutment straight Gingiva height 1.5 mm	Ti6Al4V	MU2S 1.5 GBL RC	BM6673	G
Abutment straight Gingiva height 2.5 mm	Ti6Al4V	MU2S 2.5 GBL RC	BM6674	G
Gingivaformer incl. SF MU2 Height above abutment shoulder 6 mm	Ti6Al4V	GF MU 2	BM2362	C
Localicer® incl. SF MU2 Height above abutment shoulder 6.7 mm Use with NCS Set REF 462338	Ti6Al4V	MU 2	BM2363	C
Prosthetic screw for MU2	Ti6Al4V	SFB RC	BM6176	B







ACCESSORIES FOR MULTI-UNIT ABUTMENTS

	Description	Material	Code	REF	Price cat.
	Temporary base SF MU2 sold separately	Ti6Al4V	TC MU2	BM2355	D
	Transfer straight incl. screw SFL MU2	Ti6Al4V	TS MU2	BM2356	C
	Castable for Multi-Unit incl. screw		PA MU2	BM2357	A
	Screw for TC MU2	Ti6Al4V	SF MU2	BM2358	B
	Lab analogue for Multi-Unit	Ti6Al4V	IA MU2	BM2360	B
	Hex instrument long, Ø 1.25 mm		HT 1.25	BM3022	C
	Hex instrument extralong: 45 mm, Ø 1.25 mm		HTX 1.25	BM7764	C
	Hex instrument for suprastructures, Ø 1.77 mm		HT 1.77	BM3024	C







INSTRUMENTS

	Description	Code	REF	Price cat.
	Pilot drill short/long 2.0 mm Ø	DS 2	BM1359	D
		DSL 2	BM6555	
	Pilot drill short/long 2.8 mm Ø	DS 2.8	BM1404	D
		DSL 2.8	BM6557	
	Form drill short 2.8 mm Ø	DBL 2.8	BM6033	E
	Form drill short 3.5 mm Ø	DBL 3.5	BM6034	E
	Form drill short 4.0 mm Ø	DBL 4.0	BM6035	E
	Cortical countersink 3.3	CSBL 3.3	BM6038	D
	Cortical countersink 4.1	CSBL 4.1	BM6039	D
	Cortical countersink 4.8	CSBL 4.8	BM6040	D
	Tap	TAP GBL 3.3	BM6041	D
	Tap	TAP GBL 4.1	BM6042	D
	Tap	TAP GBL 4.8	BM6043	D




GUIDE JACKET

	Description	Unit	Material	REF	Price cat.
	BFH 2.0 guide jacket 2.0 mmd	Pack of 5	Ti6Al4V	BM7100	B
	BFH 2.5 guide jacket 2.5 mmd	Pack of 5	Ti6Al4V	BM7101	B
	BFH 3.0 guide jacket 3.0 mmd	Pack of 5	Ti6Al4V	BM7102	B
	BFH 3.2 guide jacket 3.2 mmd	Pack of 5	Ti6Al4V	BM7103	B
	BFH 3.5 guide jacket 3.5 mmd	Pack of 5	Ti6Al4V	BM7104	B
	a) Length b) Height of step c) Max. Ø top d) Nominal Ø e) Ø of drilling in the drill template	5 mm 0.7 mm 3.7 - 5 mm 3 - 4.4 mm 2.05 - 3.55 mm			

ADAPTER

	Description	For	Length	Code	REF	Price cat.
	Adapter short / contra-angle	ITV GBL	22 mm	ITV S WST	BM1366	D
	Adapter long / contra-angle	ITV GBL	32 mm	ITV L WST	BM1367	D
	Adapter medium / contra-angle	ITV GBL	27 mm	ITV M WST	BM1368	D
	Ratchet adapter	Adapter for ITV GBL		IT ITV	BM1365	D
	Drill extension Contra-angle, extends by 19 mm			DX2	BM1349	D
	Universal adapter For all contra-angle instruments Use with ratchet TW2 or RAT 2, max. 30 Ncm			UAW	BM3026	E

INSTRUMENTS AND TOOLS

	Description	Type	REF	Price cat.
	Ratchet RAT 2	For all Hex instruments and insertion tools	BM1352	K
	TW2	Torque wrench, 10 - 70 Ncm. For all insertion tools, hex- and torxinstruments <i>It is recommended to have the torque ratchets recalibrated by us once a year.</i>	BM1356	S
	TT 1.25	Torx instrument (for all screws)	BM3027	C
	TT 1.25 M	Torx instrument (all screws) for contra-angle	BM3028	C
	HT 1.77	Hex instrument, long	BM3024	C
	HTX 1.77	Hex instrument, extralong	BM1070	C
	PUW1	Punch	BM3002	C

STARTER TRAY

Description	Code	REF	Price €
Adapter contra-angle short	ITV S	BM1366	
Adapter contra-angle medium	ITV M	BM1368	
Ratchet adapter for IT V	IT ITV	BM6069	
Torx instrument	TT 1.25	BM3027	
Cortical countersink 3.3	CSBL 3.3	BM6038	
Cortical countersink 4.1	CSBL 4.1	BM6039	
Cortical countersink 4.8	CSBL 4.8	BM6040	
Pilot drill	DS 2.0	BM1359	
Form drill	DBL 2.8	BM6033	
Form drill	DBL 3.5	BM6034	
Form drill	DBL 4.0	BM6035	
Tap	TAP GBL 3.3	BM6041	
Tap	TAP GBL 4.1	BM6042	
Tap	TAP GBL 4.8	BM6043	
Torque wrench	TW2	BM1356	
Starter tray empty		BM6502	upon request
Starter tray with content		SBM6008	upon request

INSTRUMENT TRAY

Description	Code	REF	Price €
Pilot drill	DS 2	BM1359	
Form drill	DBL 2.8	BM6033	
Form drill	DBL 3.5	BM6034	
Form drill	DBL 4.0	BM6035	
Standardized probe	PDG	BM1350	
Standardized probe	PDG	BM1350	
Standardized probe	PDG	BM1350	
Cortical countersink 3.3	CSBL 3.3	BM6038	
Cortical countersink 4.1	CSBL 4.1	BM6039	
Cortical countersink 4.8	CSBL 4.8	BM6040	
Tap	TAP GBL 3.3	BM6041	
Tap	TAP GBL 4.1	BM6042	
Tap	TAP GBL 4.8	BM6043	
Ratchet adapter for IT V	IT ITV	BM6069	
Adapter contra-angle short	ITV S	BM1366	
Adapter contra-angle medium	ITV M	BM1368	
Adapter contra-angle long	ITV L	BM1367	
Universal adapter	UAW	BM3026	
Punch	PUW 1	BM3002	
Torx instrument	TT 1.25	BM3027	
Drill extension	DX 2	BM1349	
Torque wrench	TW2	BM1356	
Instrument tray empty		BM6008	upon request
Instrument tray with content		SBM6008	upon request

ONEWAY BIOMED

We are certified DIN EN ISO 13485, and annex II of EEC Directive 93/42 EWG (2007).
Product dimension described in this brochure may differ from reality for technical reasons.
GBL® implants are protected by patents. GBL® is a registered trademark.
In case that implants would be reprocessed (cleaned, resterilized) infections could occur,
because no validated procedures for reprocessing are available.

(The products of this catalogue are CE marked
(class I) and CE 1936 marked (class IIa and IIb)
according to 93/42/EC Directive).

Commercial products that are not monitored
by our notified body are declared as third-party
products.

Compilation and clarification of symbols on the pack:



Batch No.



Sterilized by
gamma radiation



Non-sterile



Intended for use
by dentists or
surgeons only



Single use
product



Instruction
for use



Expiry date



Store
in a dry
place



Store tightly
keep closed



Do not use if
packing is
damaged



Do not
resterilize



Manufacturer



Production
date



Catalogue
number



Safely anti-rotational thanks to its internal precision square

Cone technology for a tight seal

Universally suitable for fixed and removable prosthodontics

The cone centers the abutment and provides 100% tightness

ONEWAY
BIOME



Dr. Ihde Dental AG
Dorfplatz 11
CH - 8737 Gommiswald / SG
Tel +41 (0)55 293 23 23
contact@implant.com
www.implant.com

Distributed by onewaybiomed GmbH
Dorfplatz 11
CH - 8737 Gommiswald / SG
Phone +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.onewaybiomed.com



Dr. Ihde Dental GmbH
Erfurter Str. 19
D - 85386 Eching / München
Tel +49 (0)89 319 761 0
info@ihde-dental.de